



MultiOx™ adjustable venturi valve

Offers the user seven oxygen concentrations



Oxygen and Aerosol Therapy ▪ Fixed Oxygen Concentration (High Flow)



MultiOx™ adjustable venturi valve

The MultiOx adjustable venturi valve offers the user seven oxygen concentrations. A large body diameter provides space for clear, widely spaced calibration and incorporates a 'pull-twist' action preventing accidental adjustment. The one-piece design eliminates the need to use more than one valve for all concentrations and there are no separate pieces to get lost or mislaid.

Ease of use

The MultiOx adjustable venturi valve is simple to use and has the following features:

- The anti tamper ring with its 'pull-twist' action is designed to prevent inadvertent adjustment
- A flow in L/min is printed next to each concentration indicating the optimum flow meter setting. These flows have been selected to ensure that the patient receives a total flow of approximately 40 L/min. This is designed to meet the patient's 'normal peak inspired flow' requirements¹

Precision

The MultiOx delivers precise oxygen concentrations in compliance with ISO Standard EN 13544-3:2001.

This accuracy is independent of the oxygen flow meter setting and is achieved using just one valve.

Comfortable for the environment

Various groups have questioned the use of PVC in medical products and it's impact on the environment. Intersurgical's aim is to reduce the environmental impact of its products and processes. The utilisation of new materials and the latest manufacturing technology has resulted in the elimination of PVC from the mask shell resulting in a reduced environmental impact.

More about
our Eco Range

www.intersurgical.com/info/eco

Incurved nose seal

Conforms to different nose shapes designed to prevent oxygen entering patient's eyes

Large, well spaced calibration settings

Easy to read, minimises errors

One piece design

No parts to change or lose

Make an enquiry



Soft face seals
Increased patient comfort

Elastic positioned under the ears
Eliminates trauma to top of ears

Pull-twist controls
Prevent inadvertent adjustment

Code	Description	Tube length	Box Qty.
1117005	MultiOx, Eco, adult, adjustable venturi mask kit with tube	1.8m	80
1118005	MultiOx, Intersurgical EcoLite, paediatric, adjustable venturi mask kit with tube	1.8m	40
1198015	Intersurgical EcoLite™, paediatric, aerosol mask		42
1188000	Eco, adult, aerosol mask		50

Reference: 1.Clinical practice in Respiratory Care. James B.Fink.



Lower environmental impact product



INTERSURGICAL
COMPLETE RESPIRATORY SYSTEMS

Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK

T: +44 (0)118 9656 300 F: +44 (0)118 9656 356 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 9001:2008, ISO 13485:2003 and ISO 14001:2015

Please think before you print

Save energy and paper.
If you must print this information sheet please print it double sided.

IS10.29 • Issue 5 04.18